

## Miami J<sup>®</sup>

Demonstrated by multiple studies to be the superior immobilizer on the market, the patented Miami J<sup>®</sup> Collar features a phenotype-driven sizing system that simplifies size selection, with 6 unique front and back components that are engineered to work optimally together for various anatomies.



### Indications for Use

- C-Spine precaution for trauma patients
- Immobilization for pre and post c-spine surgery
- Degenerative disorders
- Spinal Stenosis
- Spondylolisthesis

### Product Highlights

- A phenotype-driven sizing system — with unique fronts and backs that are bioengineered to work together for superior immobilization — makes size selection simple, while drastically reducing the margin of error in fitting a patient
- Adjustable mandibular and occipital support enables customization to disburse pressure and avoid skin breakdown
- Sorbatex™ pads that are specifically engineered to be anti-decubitus — protecting skin integrity during extended wear — and are removable for easy cleaning
- Sternal Pad disburses pressure and is removable for comfort and improved compliance when supine and/or swallowing
- Patient Assist Strap allows for flexibility and comfort without sacrificing immobilization
- Latex-free and compatible with all imaging modalities
- [Order Information](#)
- [Support Documents](#)
- [Education](#)

## Miami J Cervical Collar: Collar Sets and Accessories

Collar	Collar Set	Replacement Pads	Size
MJ-200S	MJR-200S	MJP-100	Super Short*
MJ-200L	MJR-200L	MJP-200L	Stout
MJ-250	MJR-250	MJP-250	XS
MJ-500	MJR-500	MJP-100	Large
MJ-400	MJR-400	MJP-100	Medium
MJ-300	MJR-300	MJP-100	Small

\* Super Short collars are packaged with Medium Back Panels.

## Miami J Cervical Collar: Front & Back Panels

Miami J Front	Miami J Back	Size
MJ-200SFRT	MJ-400BK	Super Short*
MJ-200LFRT	MJ-200LBK	Stout
MJ-250FRT	MJ-250BK	XS
MJ-500FRT	MJ-500BK	Large
MJ-400FRT	MJ-400BK	Medium
MJ-300FRT	MJ-300BK	Small

Miami J Front and Back Panels may be interchanged to accommodate patient anatomy, e.g., a Small Front may be used with a Large Back.

\* Super Short collars are packaged with Medium Back Panels.

### Please note:

Choose collar size based on patient neck height, not patient body size.

Each collar has a bioengineered front with a specific back designed to work together for superior immobilization.